

EMPTY BAGS

MEDICAL AND PHARMACEUTICAL
SOLUTIONS

FLEXIBILITY AND CUSTOMIZATION
CONTRACT MANUFACTURING EXPERTISE

HAEMOPHARM
RESEARCH AND INNOVATION

PVC-FREE AND PVC

MEDICAL AND PHARMACEUTICAL PACKAGING

QUALITY FLEXIBILITY AND CUSTOMIZATION

Study, design and manufacture of empty bags
for medical and pharmaceutical solutions

QUALITY IS NEVER A COMPROMISE



PVC-FREE BAGS


VOLUMES

single and multi-compartment bags
from ml 50 to ml 6000
other volumes are available on request

CHARACTERISTICS

from 1 to 5 tubes having each
length and diameter as desired
hot-printed label one
or two colours

PVC BAGS



USE
i.v. and dialysis solutions
antibiotics
medical device fluids
clinical nutrition
cryopreservation

AUTOMATIC MACHINE FOR MEDICAL-GRADE BAGS



The tomorrow's technology
available today for our
customers

HAEMOPHARM HEALTHCARE
I - 20135 MILANO - V.le CIRENE, 1
Tel. +39 02 5513513 - Fax +39 02 55191982
info@haemopharm.it
www.haemopharm.it