

The background of the entire page is a light blue gradient with numerous water droplets of various sizes, some in sharp focus and others blurred, creating a sense of freshness and purity.

**SANO<sub>3</sub>**

**KIT OF  
OZONE THERAPY  
FOR MAJOR INFUSION  
AUTOHEMO THERAPY PROTOCOL**

**HAEMOPHARM**  
RESEARCH AND INNOVATION



### **CLOSED CIRCUIT**

SanO3 is a closed circuit system that allows to carry out all the procedures of the Major Autohemo Infusion Therapy in conditions of extremely high safety, avoiding risks of environmental contamination.

### **NO NEEDLES**

SanO3 system does not require other needles except the sample one.

All the procedures of the anticoagulant and O<sub>2</sub> O<sub>3</sub> filling, are taking place with no needle syringes, so without risks to contract infections for accidental pricks.

### **PHTHALATES FREE MATERIALS**

All the components in SanO3 system are certified by the major International Pharmacopoeia and subjected to a selection taking into consideration the highest qualitative requirements. The entire system is phthalates free.

### **AUTOMATIC CONNECTIONS**

SanO3 system is provided with vanguard's advanced connections with automatic open-closed ports, assuring hermetic closing before and after use.

### **ABSENCE OF BLOOD CLOTS AND GAS EMBOLUS**

SanO3 permits to carry out all the necessary procedures for the Major Autohemo Infusion Therapy treatment without risks of blood clots and gas embolus.

**FROM HAEMOPHARM'S RESEARCH IS TURNED OUT  
A NEW SPECIFIC KIT FOR THE INTRAVENOUS MAJOR AUTOHEMO  
INFUSION THERAPY, ABSOLUTELY INNOVATIVE AND SAFE.**

## THE KIT



## TECHNICAL CHARACTERISTICS

A 400 ml. Bag in PVC free from phthalates  
Containing 35 ml. of anticoagulant ACD-A solution  
Complete of automatic connections and set for blood sample;  
Sterile;

Blood transfusion set with 175 micron filter  
Complete of Y injection port and Luer Lock connection;  
Sterile;

PATENT PENDING

CE0373

*All the information contained in this brochure  
and the described product can be changed in any  
moment by Haemopharm Healthcare without notice.*

**SANO<sub>3</sub>**

Manufactured by:

**PAOLO GOBBI FRATTINI**  
I - 23030 Tovo di S. Agata (SO)  
Str. provinciale per Vervio, 8  
Tel. +39 0342 771017 - Fax +39 0342 771018

Distributed by:

**HAEMOPHARM HEALTHCARE - MILANO**  
[info@haemopharm.it](mailto:info@haemopharm.it)  
[www.haemopharm.it](http://www.haemopharm.it)